

# EXO GAG

## SALIVA

Next generation of Exosome isolation  
for biomarker discovery in liquid biopsy

### EXOSOME PURIFICATION KIT FOR SALIVA

#### USER GUIDE

#### STORAGE

---

All components can be stored at room temperature.

#### PRODUCT COMPONENT

---

EXO GAG - Saliva™ isolation kit 20ml (20-10 samples). 1 x User guide.  
Not supplied: 1.5ml microcentrifuge collection tubes.

**EXO GAG - SALIVA is a specific, quick and inexpensive method to optimize the process of exosome isolation.**

#### PRODUCT INFORMATION

---

The EXO GAG - Saliva precipitation reaction is based on the interaction between the precipitation solution and glycosaminoglycans exosomes surface. After a simple incubation, the exosomes can be isolated by a short centrifugation.

#### WHAT IS EXO GAG - SALIVA

---

EXO GAG - Saliva technology allows the isolation of exosomes from saliva based on the affinity of exosomes to the precipitation reagent, which allows their isolation from a complex sample. EXO GAG - Saliva is a patented exosome purification method that allows the isolation of exosomes from a saliva sample with a minimal amount of co-precipitated material, such as protein or genetic material (DNA, RNA or microRNA) which makes it an ideal product for biomarkers research and their transfer to the clinic.



**High sensibility and specificity**


**Fast**

**Inexpensive**

**Small sample needed**

**No specific equipment needed**

SAMPLE	SAMPLE, vol.	EXOAG - SALIVA, vol.
Saliva	500µl	1ml
Saliva	1ml	2ml



## PROTOCOL

### EXOAG - Saliva Exosomes Isolation Protocol.

- 1. Collect saliva sample.** Samples can be frozen until the moment they are used; if the samples have been frozen, thaw and temper them before processing.
- 2. Centrifuge the sample** at 2000 x g for 5' to remove cells and cell debris.
- 3. Transfer the supernatant** to a new tube and discard the pellet of possible cell debris.
- 4. Add the volume of sample** to isolate exosomes to a new tube and add twice the volume of EXOGAG - Saliva precipitation reagent, as shown in the table.
- 5. Mix the sample** and EXOGAG - Saliva precipitation reagent by inverting the tube or vortexing to homogenize the final solution (the solution will have a characteristic blue colour).
- 6. Incubate the sample** for 5' at 4°C.
- 7. Centrifuge the sample** at 3000 x g; 15' at 4°C.
- 8. Remove the supernatant** being careful not to remove the pellet containing the exosomes (this pellet will be dark blue).
- 9. Resuspend the exosomes** in the appropriate buffer (repeatedly pipetting up and down), depending on the technique.